

## N O T I C E

THIS DOCUMENT HAS BEEN REPRODUCED FROM  
MICROFICHE. ALTHOUGH IT IS RECOGNIZED THAT  
CERTAIN PORTIONS ARE ILLEGIBLE, IT IS BEING RELEASED  
IN THE INTEREST OF MAKING AVAILABLE AS MUCH  
INFORMATION AS POSSIBLE



National Aeronautics and  
Space Administration

CR 160895

c.4

**Lyndon B. Johnson Space Center**  
Houston, Texas 77058

# ELECTRON MICROPROBE MINERAL ANALYSIS GUIDE

Roy W. Brown  
Lockheed Engineering and  
Management Services Co., Inc.  
Planetary and Earth Sciences Dept.  
Lyndon B. Johnson Space Center  
Houston, Texas 77058

NAS9-15800

September 1980



Prepared for  
LYNDON B. JOHNSON SPACE CENTER  
Houston, Texas 77058

(NASA-CR-160895) ELECTRON MICROPROBE  
MINERAL ANALYSIS GUIDE (Lockheed Engineering  
and Management) 187 p HC A09/MF A01

N81-14018

CSCL 07D

Unclass

G3/25

29504



1. Report No. NASA CR-	2. Government Accession No. —	3. Recipient's Catalog No.	
4. Title and Subtitle  ELECTRON MICROPROBE MINERAL ANALYSIS GUIDE		5. Report Date September 1980	
		6. Performing Organization Code —	
7. Author(s) Roy W. Brown		8. Performing Organization Report No.	
9. Performing Organization Name and Address Lockheed Engineering and Management Services Co., Inc. Planetary and Earth Sciences Dept. Lyndon B. Johnson Space Center, Houston, TX		10. Work Unit No.	
		11. Contract or Grant No. NAS9-15800	
12. Sponsoring Agency Name and Address Lyndon B. Johnson Space Center Houston, Texas 77058 Dr. M. B. Duke - Technical Monitor		13. Type of Report and Period Covered	
		14. Sponsoring Agency Code —	
15. Supplementary Notes  Presented as a laboratory guide to X-ray mineral analysis using the electron microprobe.			
16. Abstract  Electron Microprobe Mineral Analysis Guide is a compilation of X-ray tables and spectra recorded from various mineral matrices. Spectra were obtained using an Applied Research Laboratories (A.R.L.) electron microprobe, model EMX-SM equipped with LiF geared, curved crystal X-ray spectrometers, utilizing typical analytical operating conditions: 15 Kv acceleration potential, 0.02 micro-Ampere sample current as measured on a clinopyroxene standard (CP19). Tables and spectra are presented for the majority of elements, fluorine through uranium, occurring in mineral samples from lunar, meteoritic and terrestrial sources. Tables for each element contain relevant analytical information, i.e., analyzing crystal, X-ray peak, background and relative intensity information, X-ray interferences and a section containing notes on the measurement. Originally intended to cover silicates and oxide minerals the tables and spectra have been expanded to cover other mineral phases. Electron Microprobe Mineral Analysis Guide is intended as a "spectral base" to which additional spectra can be added as the analyst encounters new mineral matrices.			
17. Key Words (Suggested by Author(s))  X-ray, electron microprobe, quantitative mineral analysis		18. Distribution Statement	
19. Security Classif. (of this report) N/A	20. Security Classif. (of this page) N/A	21. No. of Pages 184	22. Price*

\*For sale by the National Technical Information Service, Springfield, Virginia 22161

## PREFACE

The object of this report is to provide a laboratory guide of X-ray spectra used in quantitative electron microprobe analysis of minerals. This compilation of X-ray spectra covers the majority of elements encountered during mineral analysis from the element fluorine (atomic number 9) through uranium (atomic number 92). Spectra were obtained from lunar, meteoritic and terrestrial mineral samples using an A.R.L. EMX-SM electron microprobe. Hopefully, this report would provide in a single source, a reasonable assessment of interference problems which confront the microprobe analyst.

PRECEDING PAGE BLANK NOT FILMED

# LIST OF TABLES

TABLE NO. 9	FLUORINE, F	TABLE NO. 66	DYSPROSIUM, Dy
TABLE NO. 11	SODIUM, Na	TABLE NO. 67	HOLMIUM, Ho
TABLE NO. 12	MAGNESIUM, Mg	TABLE NO. 68	ERBIUM, Er
TABLE NO. 13	ALUMINUM, Al	TABLE NO. 69	THULIUM, Tm
TABLE NO. 14	SILICON, Si	TABLE NO. 70	YTTERBIUM, Yb
TABLE NO. 15	PHOSPHORUS, P	TABLE NO. 71	LUTETIUM, Lu
TABLE NO. 16	SULFUR, S	TABLE NO. 78	PLATINUM, Pt
TABLE NO. 17	CHLORINE, Cl	TABLE NO. 90	THORIUM, Th
TABLE NO. 19	POTASSIUM, K	TABLE NO. 92	URANIUM, U
TABLE NO. 20	CALCIUM, Ca	TABLE 001	X-RAY OVERLAP AMONG TRANSITION ELEMENTS
TABLE NO. 22	TITANIUM, Ti		
TABLE NO. 23	VANADIUM, V	TABLE 002	X-RAY OVERLAP AMONG RARE EARTH ELEMENTS
TABLE NO. 24	CHROMIUM, Cr		
TABLE NO. 25	MANGANESE, Mn		
TABLE NO. 26	IRON, Fe		
TABLE NO. 27	COBALT, Co		
TABLE NO. 28	NICKEL, Ni		
TABLE NO. 38	STRONTIUM, Sr		
TABLE NO. 40	ZIRCONIUM, Zr		
TABLE NO. 56	BARIUM, Ba		
TABLE NO. 57	LANTHANUM, La		
TABLE NO. 58	CERIUM, Ce		
TABLE NO. 59	PRASEODYMIUM, Pr		
TABLE NO. 60	NEODYMIUM, Nd		
TABLE NO. 62	SAMARIUM, Sm		
TABLE NO. 63	EUROPIUM, Eu		
TABLE NO. 64	GADOLINIUM, Gd		
TABLE NO. 65	TERBIUM, Tb		

## Introduction

Electron Microprobe Mineral Analysis Guide (EMMAG) is a compilation of X-ray tables and spectra recorded from various mineral matrices. The Guide provides for each element complete documentation of the spectrum in the area of the selected peak and records optimum positions in the X-ray continuum for the measurement of background, and lists interferences, with both peak and background measurements, that may be encountered in mineral analysis.

Initially EMMAG was directed toward silicate and/or oxide mineral analysis and the X-ray lines listed in each table reflect both this and the emphasis of our work which at that time was primarily on lunar samples. To make the guide more general, other mineral matrices have been added as analytical problems arose. Deciding which X-ray lines to include is difficult, at times subjective, and very much dependent on my knowledge of mineral matrices and on the accuracy of the intensity information listed in "X-ray Emission and Absorption Edge Wavelengths and Interchange Settings for LiF Geared Curved Crystal Spectrometer" by E.W. White and G.G. Johnson, Jr., 2nd edition. All of the X-ray line information, i.e., wavelength, KeV, order of reflection (N), intensity (I), etc. is derived from this source. There are bound to be errors of omission in striving to list only the X-ray lines one is most likely to encounter. EMMAG is a guide to X-ray interference and background problems. In laboratory usage it is best kept in a looseleaf notebook and supplemented as one encounters new problems and matrices during the course of microprobe analysis. EMMAG is in part a condensation of

X-ray wavelength information contained in the White and Johnson, Jr. listing: between the  $\pm\Delta$  Bkg for  $\text{Si}_{\text{K}\alpha 1,2}$  are listed 156 x-ray lines - EMMAG lists only 9. The 9 X-ray lines selected are believed to be those most likely to cause problems during mineral analysis.

### Tables

X-ray tables are listed by atomic number, i.e., information for fluorine is listed in Table 9. In this manner other elements can be added without disrupting the order of the guide. Each table is divided into four sections. The first section gives the X-ray line information for the element: analytical line ( $\text{K}\alpha 1,2$ ,  $\text{L}\alpha 1$ ,  $\text{M}\alpha$ ), analyzing crystal (LiF, ADP, PET, RbAP), spectrometer setting for a given crystal, and background setting used in this laboratory. The next section, Elemental Scans, lists the minerals for which spectra have been recorded for the particular X-ray line, the wt. % of the element or oxide in the mineral and the background setting which may vary depending on the mineral matrix. The third section, Interferences, lists X-ray lines which occur within the  $\pm\Delta$  Bkg limits of the analysed line, the reflection order (N), the intensity relative to the strongest line (the unresolved  $\text{K}\alpha 2$ ), the spectrometer odometer setting for the analyzing crystal, and the energy (KeV) of the X-ray line. The spectrometers are the LiF geared type, i.e., with a LiF crystal in use the odometer reads directly in angstroms. Other analyzing crystals indicate an interchange value. The fourth section, Notes, contains comments pertaining to the interferences which occur and how to deal with them plus additional information pertinent to the analytical problem. Typically the use of a Pulse Height Analyzer (PHA) is discussed, the PHA being the main form of electronic pulse discrimination used with proportional X-ray counters.

### Spectral Scans

The X-ray scans are an important part of the guide because they give relative intensities, position, and crystal line broadening effects of X-ray lines derived from various mineral matrices. Each scan was recorded on an A.R.L. EMX-SM electron microprobe in the Geochemistry Branch of the Planetary and Earth Sciences Division, Johnson Space Center. The microprobe has a  $52.5^\circ$  X-ray take-off angle and was operated with an acceleration potential of 15 KV and sample current of 0.02 microamperes measured on a clinopyroxene standard (CP19). Crystal spectrometers have a receiver slit width of 0.015" as measured with a feeler gage. Each scan was obtained by motor driving the crystal spectrometer at a speed of  $0.02^\circ/\text{min}$ . synchronously with a Texas Instruments Servo/riter II recorder driven at 1"/min. Recorder output is then  $0.02^\circ/\text{inch}$  of chart paper. Full scale intensity is either 300 or 100 cps as labeled. Due to low values for full scale intensity major X-ray peaks are typically off scale and therefore the associated peak position is not accurately determined. More important than accurate peak position is the relative position of the minor X-ray peaks, the extent of major peak tailing and the slope of the X-ray continuum which together provide information on the interferences present as well as optimum background positions. Each scan is labeled with element symbol, mineral matrix scanned, analyzing crystal utilized, counts per second (cps) full scale deflection and wt. % of element or oxide present in the sample.

### Historical Note

Originally X-ray tables and spectra were produced for our laboratory use only; later on I thought they might be of interest to a wider audience. Because the initial intent of production was limited (1 copy) and the data collection time period (approximately nine years), some inconsistencies were propagated. Most notable among these is the use of subscripts. Therefore, let the reader be advised that  $K_{\alpha 1,2}$  and  $K_{\alpha 1,2}$  are the same X-ray line. Other inconsistencies and/or outright errors which may exist are not intentional and are solely the responsibility of the author.

### Acknowledgments

I would like to thank those NASA staff scientists and visiting scientists to the Johnson Space Center who encouraged me to formalize the guide book and have it printed for circulation. I would also like to thank "LEMSCO", Lockheed Engineering and Management Services Co., Inc., for granting me the time to work on the project and Lina Romero for doing all the typing of tables and labels without complaint.

FLUORINE



TABLE 9


ELEMENT: F  
 ANALYTICAL LINE: K $\alpha$   
 CRYSTAL: RbAP  
 SPECTROMETER SETTING: 2.8241  
 BACKGROUND SETTING:  $\pm 0.063$   
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT. %</u>	<u>BACKGROUND SETTING</u>
Topaz	~19.0	$\pm 0.063$

## INTERFERENCES:

<u>LINE</u>	<u>N</u>	<u>I</u>	<u>RbAP (100)</u>	<u>KeV</u>
P SK $\alpha$ 6	3	.01	2.8040	2.04
F SK $\alpha$ ''	1	30	2.8056	0.68
F SK $\alpha$ '	1	35	2.8086	0.68
P SK $\alpha$ 5	3	.01	2.8095	2.04
Co L1	1	9	2.8197	0.68
F K $\alpha$	1	100	2.8241	0.68
P SK $\alpha$ 4	3	5	2.8251	2.03
P SK $\alpha$ 3	3	5	2.8287	2.03
Ce MZ	1	.01	2.8287	0.68
P SK $\alpha$ '	3	.5	2.8353	2.02
P K $\alpha$ 1,2	3	150	2.8478	2.01
Ca K $\beta$ 1,3	6	15	2.8577	4.01
Dy M $\beta$	2	45	2.8848	1.33

NOTES: Fluorine is typically of interest in apatites, mica and topaz. Interferences are P and Ca, both major constituents of apatites. Both are higher order X-ray lines and a PHA can be used effectively to discriminate against them. From experience with topaz analysis, it is extremely important to match the standard and sample matrix as closely as possible due to uncertainties in the absorption correction.



95

ORIGINAL PAGE IS  
OF POOR QUALITY

SODIUM

TABLE 11

ELEMENT: Na

ANALYTICAL LINE K $\alpha$ 1,2

CRYSTAL: RbAP

SPECTROMETER SETTING: 1.8360

BACKGROUND SETTING:  $\pm 0.062$

ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT. % OXIDE</u>	<u>BACKGROUND SETTING</u>
Oligoclase	8.68	$\pm 0.062$
Kaersutite	2.80	$\pm 0.062$
Apatite Zirconium Zinc	0.30	

## INTERFERENCE:


	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>RbAP (100)</u>	<u>KeV</u>
P	KB	2	3	1.7869	2.139
P	KB1	2	3	1.7893	2.136
Fe	K $\alpha$ 1,2	6	150	1.7919	6.398
Zr	L $\beta$ 1	2	45	1.7993	2.124
Na	SK $\alpha$ 6	1	1	1.8014	1.061
Na	SK $\alpha$ 5	1	1	1.8067	1.058
Na	SK $\alpha$ 4	1	10	1.8168	1.052
Na	SK $\alpha$ 3	1	10	1.8198	1.050
Na	SK $\alpha$ 1	1	3	1.8247	1.047
Na	K $\alpha$ 1,2	1	100	1.8360	1.041
Zn	L $\beta$ 1	1	26	1.8472	1.034

<u>LINE</u>	<u>N</u>	<u>I</u>	<u>RbAP (100)</u>	<u>KeV</u>
Ba L $\beta$ 2,15	5	20	1.8532	5.156
Cu L $\beta$ 3,4	1	1	1.8686	1.023
Zr L $\alpha$ 1	2	100	1.8716	2.042
P SK $\alpha$ 4	2	5	1.8834	2.029
P SK $\alpha$ 3	2	5	1.8858	2.027
Zn L $\alpha$ 1,2	1	100	1.8890	1.012
P K $\alpha$ 1,2	2	150	1.8985	2.013
Ca KB1,3	4	15	1.9051	4.012

NOTES: Minerals with high P concentrations might cause high background counts.

Interference with Na<sub>K $\alpha$ 1,2</sub> peak measurements could occur in a Zn matrix (ZnL $\beta$ 1). Energies of the ZnL $\beta$ 1 (1.03 Kev) and Na<sub>K $\alpha$ 1,2</sub> (1.04 Kev) peaks are essentially identical and not separable by PHA methods.


Background measurements can be difficult in a high P matrix (PK $\alpha$ 1,2 2nd order line) but P line should be excluded with use of PHA window. A very high Fe concentration matrix could conceivably cause a problem but could also be discriminated against with use of a PHA window. Zr has two X-ray peaks (ZrL $\alpha$ 1, ZrL $\beta$ 1) in this region — both should be excluded by PHA methods. Au M $\alpha$ 1,2 also occurs in this region but is commonly not a problem.




13

FEET  
LEFT

Na kaersutite  
RDAP 300 cps.  
Na<sub>2</sub>O 2.80 wt. %  
FeO 11.87 wt. %





ORIGINAL PAGE IS  
OF POOR QUALITY





57





MAGNESIUM

TABLE 12

ELEMENT: Mg

ANALYTICAL LINE:  $K\alpha_{1,2}$

CRYSTAL: RbAP

SPECTROMETER SETTING: 1.5246

BACKGROUND SETTING:  $\pm 0.093$

ELEMENTAL SCANS:


<u>MINERAL</u>	<u>WT. % OXIDE</u>	<u>BACKGROUND SETTING</u>
Olivine, Marjalahti	47.42	$\pm 0.093$
Lake Co. plagioclase	$\sim 600$ ppm	
Kaersutite	11.36	$\pm 0.093$
Ilmenite	9.77	$\pm 0.093$
Chromite	8.09	$\pm 0.093$

## INTERFERENCES:

<u>LINE</u>	<u>N</u>	<u>I</u>	<u>RbAP (100)</u>	<u>KeV</u>
Cr $K\alpha_{1,2}$	4	150	1.4126	5.411
Ca $K\beta_{1,3}$	3	15	1.4288	4.012
Cl $K\alpha_{1,2}$	2	150	1.4579	2.621
Mn $K\beta_{1,3}$	5	20	1.4723	6.489
Ba $L\beta_{2,15}$	4	20	1.4825	5.156
Fe $K\alpha_{1,2}$	5	150	1.4932	6.398
Co $K\beta_{1,3}$	6	20	1.4991	7.648

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>RbAP</u> <u>(100)</u>	<u>KeV</u>
Mg	SK $\alpha$ 4	1	8	1.5119	1.264
Mg	SK $\alpha$ 3	1	8	1.5144	1.262
Mg	SK $\alpha$ 1	1	2	1.5181	1.259
Mg	K $\alpha$ 1,2	1	150	1.5246	1.253
Ni	K $\alpha$ 1,2	6	150	1.5346	7.471
Ca	SK $\alpha$ 4	3	2	1.5431	3.715
V	K $\alpha$ 1,2	4	150	1.5445	4.949
Ca	SK $\alpha$ 3	3	2	1.5447	3.711
S	K $\beta$ <sub>x</sub>	2	1	1.5487	2.468
Ti	K $\beta$ 1,3	4	20	1.5501	4.931
S	K $\beta$ 1	2	7	1.5513	2.464
Ba	L $\beta$ 3	4	6	1.5516	4.926
Ca	K $\alpha$ 1,2	3	150	1.5536	3.690
Ba	L $\beta$ 4	4	5	1.5756	4.851
Ba	L $\beta$ 1	4	5	1.5836	4.827
K	K $\beta$ 1,3	3	15	1.5973	3.589
Cr	K $\beta$ 1,3	5	18	1.6069	5.946


NOTES: Problems with Mg measurement might occur with a high Ca, Ni or V matrix. All 3 lines are on the high angle side of the Mg<sub>K $\alpha$ 1,2</sub> peak. Ni<sub>K $\alpha$ 1,2</sub> 6th and V<sub>K $\alpha$ 1,2</sub> 4th are closest to the Mg peak but do not occur in sufficient concentration in a typical silicate matrix to be a problem, however, their energies are sufficiently different from the Mg line so that they could be discriminated against with a PHA. A more common problem is a high Ca matrix (plagioclase). Fortunately the most intense Ca line (Ca<sub>K $\alpha$ 1,2</sub> 3rd) can also be removed with a PHA. (See scan labeled Lake Co. Plagioclase.) In choosing a background delta value, consideration should be given to the presence of numerous possible X-ray lines occurring in the Mg<sub>K $\alpha$ 1,2</sub> region.







Mg Lake Co. Plagioclase  
RbAp 300 cps.  
MgO ~0.09 wt. %


E = 1.30  
 $\Delta E = 1.90$











Mg chromite  
RBP 300 cps.  
Mg 8.09 wt. %

94  
FEET  
LEFT

ALUMINUM

TABLE 13

ELEMENT: Al

ANALYTICAL LINE:  $K\alpha_{1,2}$

CRYSTAL: RbAP; ADP

SPECTROMETER SETTING: 1.2856 (RbAP); 3.1562 (ADP)

BACKGROUND SETTING:  $\pm 0.074$  (RbAP);  $\pm 0.060$  (ADP)

ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT. % OXIDE</u>	<u>BACKGROUND SETTING</u>	
		(RbAP)	(ADP)
Bytownite	32.05	$\pm 0.074$	$\pm 0.060$
Kaersutite	14.89	$\pm 0.074$	$\pm 0.060$
Chromite	12.77	$\pm 0.074$	$\pm 0.060$
Al-Fe Spinel	10.00	$\pm 0.074$	
Ilmenite	0.26		$\pm 0.010$
Armstrongite	$\approx 1.60$		$\pm 0.010$


## INTERFERENCES:

<u>LINE</u>	<u>N</u>	<u>I</u>	<u>RbAP</u>	<u>ADP</u>	<u>KeV</u>
V $K\alpha_{1,2}$	3	150	1.1584	2.8438	4.949
Ti $K\beta_{1,3}$	3	20	1.1626	2.8541	4.931
Ba $L\beta_3$	3	6	1.1637	2.8569	4.926
Mn $K\beta_{1,3}$	4	20	1.1778	2.8916	6.489
Ba $L\beta_1$	3	50	1.1877	2.9157	4.827


	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>RbAP</u>	<u>ADP</u>	<u>KeV</u>
	Cu K $\alpha$ 1,2	5	150	1.1884	2.9175	8.040
	Fe K $\alpha$ 1,2	4	150	1.1946	2.9327	6.398
	La L $\alpha$ 1	3	100	1.2328	3.0264	4.650
	La L $\alpha$ 2	3	10	1.2372	3.0374	4.633
	Co K $\beta$ 1,3	5	20	1.2492	3.0669	7.648
	Ti K $\alpha$ 1,2	3	150	1.2716	3.1218	4.508
Al	SK $\alpha$ 4	1	4	1.2749	3.1300	1.499
Al	SK $\alpha$ 3	1	8	1.2774	3.1360	1.496
	Ni K $\alpha$ 1,2	5	150	1.2788	3.1395	7.471
Al	SK $\alpha$ 1	1	2	1.2802	3.1428	1.493
	Ba L $\alpha$ 1	3	100	1.2838	3.1516	4.465
	Cr K $\beta$ 1,3	4	18	1.2855	3.1560	5.946
Al	K $\alpha$ 1,2	1	150	1.2856	3.1562	1.486
	Cu K $\beta$ 1,3	6	20	1.2877	3.1612	8.904
	Cu K $\beta$ 3	6	6	1.2880	3.1621	8.901
	Ba L $\alpha$ 2	3	10	1.2882	3.1625	4.450
	Mn K $\alpha$ 1,2	4	150	1.2968	3.1836	5.894
	Zn K $\alpha$ 1,2	6	150	1.3286	3.2616	8.630
	Fe K $\beta$ 1,3	5	20	1.3539	3.3239	7.057
	Co K $\alpha$ 1,2	5	150	1.3799	3.3875	6.924
	Ba L $\alpha$ 1	4	5	1.3821	3.3931	5.53
	Ni K $\beta$ 1,3	6	20	1.3875	3.4063	8.263
	V K $\beta$ 1,3	4	20	1.4086	3.4580	5.426
	Cr K $\alpha$ 1,2	4	150	1.4126	3.4680	5.411

NOTES: Measurement of  $Al_{K\alpha 1,2}$  can be difficult due to the numerous peaks occurring in the region of the Al peak. Only one X-ray line ( $Cr_{K\beta 1,3^{4th}}$ ) interferes directly with the peak measurement and it can be excluded with a PHA window.

$Ti_{K\alpha 1,2^{3rd}}$  can also interfere with an Al determination, due to its close proximity to the Al peak (see ilmenite scan). Fortunately a PHA window can also be used to minimize the Ti peak. Other elements occur in sufficiently low concentration in typical silicate matrices and also are high order lines and therefore less likely to present a problem.









ORIGINAL PAGE IS  
OF POOR QUALITY









98

Al kaersutite  
ADP 300 cps.  
Al<sub>2</sub>O<sub>3</sub> 14.69 wt. %







ORIGINAL PAGE IS  
OF POOR QUALITY



Al chromite  
ADP 300 cps.  
Al<sub>2</sub>O<sub>3</sub> 12.77 wt. %






27

FEET  
LEFT

h

Al armalcolite

ADP 300 cps.

Al<sub>2</sub>O<sub>3</sub> 1.28 wt. %

SILICON



TABLE 14

ELEMENT: Si

ANALYTICAL LINE: K $\alpha$ 1,2

CRYSTAL: ADP

SPECTROMETER SETTING: 2.6969

BACKGROUND SETTING:  $\pm 0.060$

ELEMENTAL SCANS:

<u>LINE</u>	<u>WT. % OXIDE</u>	<u>BACKGROUND SETTING</u>
CP19	47.44	$\pm 0.060$

Fe

## INTERFERENCES:


	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>ADP (101)</u>	<u>KeV</u>
Fe	K $\beta$ 1,3	4	20	2.6591	7.057
Y	Ln	1	1	2.6644	1.761
Si	SK $\alpha$ 4	1	3	2.6743	1.754
Si	SK $\alpha$ 3	1	6	2.6781	1.752
Si	K $\alpha$ 1,2	1	150	2.6969	1.739
Co	K $\alpha$ 1,2	4	150	2.7100	6.924
Zn	K $\alpha$ 1,2	5	150	2.7180	8.630
Ba	L $\beta$ 2,15	3	20	2.7297	5.156
La	L $\beta$ 3	3	6	2.7367	5.143


NOTE: No interference with peak measurement in typical silicate matrix. Avoid Fe<sub>K $\beta$ 1,3</sub> peak when choosing background positions.

19

55

S1 clinopyroxene (CP19)  
ADP 300 cps.  
SiO<sub>2</sub> 48.39 wt. %





Si clinopyroxene (CP19)  
PET 300 cps.  
SiO<sub>2</sub> 48.39 wt. %




TABLE 19

ELEMENT: K

ANALYTICAL LINE:  $K\alpha_{1,2}$ 

CRYSTAL: ADP, PET

SPECTROMETER SETTING: 1.4163(ADP); 1.7238(PET)

BACKGROUND SETTING:  $\pm 0.070$ 

ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT. % OXIDE</u>	<u>BACKGROUND SETTING</u>
Orthoclase	15.36	$\pm 0.070$
Kaersutite	1.59	$\pm 0.070$

INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>ADP</u> <u>(101)</u>	<u>PET</u> <u>(002)</u>	<u>KeV</u>
	Co $K\alpha_{1,2}$	2	150	1.3550	1.6492	6.92
K	$SK\alpha_4$	1	3	1.4065	1.7118	3.34
K	$SK\alpha_3$	1	3	1.4080	1.7137	3.33
K	$K\alpha_{1,2}$	1	150	1.4163	1.7238	3.31
	Mn $K\beta_5$	2	.03	1.4359	1.7476	6.53
	Mn $K\beta_{1,3}$	2	20	1.4458	1.7597	6.49
	Fe $K\alpha_{1,2}$	2	150	1.4663	1.7847	6.40

NOTES: By using a background setting of  $\pm 0.070$  the possible background interference from  $Mn_{K\beta_{1,3}}^{2nd}$  and  $Fe_{K\alpha_{1,2}}^{2nd}$  should be avoided.

(See scans.)

PHOSPHORUS

TABLE 15

ELEMENT P

ANALYTICAL LINE:  $K\alpha_{1,2}$

CRYSTAL: PET

SPECTROMETER SETTING: 2.8364

BACKGROUND SETTING:  $\pm 0.060$

ELEMENTAL SCANS:


<u>MINERAL</u>	<u>WT. % OXIDE</u>	<u>BACKGROUND SETTING</u>
Apatite	41.29 $P_2O_5$	$\pm 0.060$
Schreibersite	$\sim 15.5$ % P	$\pm 0.060$
Diopside		

## INTERFERENCES:

	<u>N</u>	<u>I</u>	<u>ADP</u> <u>(101)</u>	<u>PET</u> <u>(002)</u>	<u>KeV</u>
Ni $K\beta_{1,3}$	4	20	2.2708	2.7639	8.26
Pt $M\alpha_1$	1	100	2.2884	2.7853	2.05
P $SK\alpha_6$	1	.01	2.2946	2.7927	2.04
Zr $L\alpha_1$	1	100	2.2973	2.7961	2.04
P $SK\alpha_5$	1	.01	2.2991	2.7983	2.04
Zr $L\alpha_2$	1	10	2.3000	2.7995	2.04
P $SK\alpha_4$	1	5	2.3118	2.8137	2.03
P $SK\alpha_3$	1	5	2.3148	2.8173	2.03
P $SK\alpha_1$	1	.5	2.3202	2.8239	2.02
Ca $K\beta_5$	2	.01	2.3271	2.8323	4.03
P $K\alpha_{1,2}$	1	150	2.3304	2.8364	2.01
Cu $K\alpha_{1,2}$	4	150	2.3340	2.8407	8.04
Ca $K\beta_{1,3}$	2	15	2.3385	2.8463	4.01
Ge $K\alpha_{1,2}$	5	150	2.3754	2.8911	9.87
Mn $K\alpha_{1,2}$	3	150	2.3877	2.9061	5.89


NOTES: The  $Ca_{K\beta_{1,3}}$  2nd interferes with the  $P_{K\alpha_{1,2}}$  measurement, however, it does not usually present a problem in silicate analysis because the Ca line is a 2nd order line and requires  $\sim 20.0$  wt. % CaO concentration before the line is observed (see diopside scan). Therefore, Ca would typically only be a problem in minerals such as apatite or a high Ca glass. Where it is a problem it can be discriminated against with a properly set PHA window.






ORIGINAL PAGE IS  
OF POOR QUALITY





P diopside  
PET 100 cps.  
CaO 24.7 wt. %



SULFUR

TABLE 16

ELEMENT : S

ANALYTICAL LINE: K $\alpha$ 1,2

CRYSTAL: ADP, PET

SPECTROMETER SETTING: 2.0334(ADP); 2.4749(PET)

BACKGROUND SETTING:  $\pm 0.065$ 

ELEMENTAL SCANS:

<u>MINERAL</u>	<u>wt. % S</u>	<u>BACKGROUND SETTING</u>
Troilite	36.47	$\pm 0.065$
Barite	13.74	$\pm 0.040$
Co		

INTERFERENCES :


<u>LINE</u>	<u>N</u>	<u>I</u>	<u>ADP (101)</u>	<u>PET (002)</u>	<u>KeV</u>
Fe K $\beta$ 5	3	.03	1.9802	2.4102	7.11
Nb L $\beta$ 2,15	1	1	1.9822	2.4126	2.37
Hg L $\beta$ 1	5	50	1.9843	2.4151	11.82
Tl M $\beta$	1	55	1.9864	2.4177	2.36
Pt L $\alpha$	4	100	1.9876	2.4192	9.44
Sr K $\alpha$ 1,2	6	150	1.9906	2.4227	14.14
Bi L1	4	3	1.9923	2.4248	9.42
Fe K $\beta$ 1,3	3	20	1.9943	2.4273	7.05
Sb L $\delta$ 4	2	.1	1.9980	2.4318	4.70
Pb M $\alpha$ 1	1	100	2.0004	2.4348	2.35
As K $\beta$ 1	5	15	2.0004	2.4350	11.72
As K $\beta$ 3	5	7	2.0016	2.4362	11.72
Pt L $\alpha$ 2	4	10	2.0047	2.4399	9.36
Pb M $\alpha$ 2	1	100	2.0054	2.4407	2.34
W L $\delta$ 3	5	2	2.0095	2.4458	11.67
Nb L $\beta$ 3	1	3	2.0096	2.4459	2.33
Pt M $\delta$	1	3	2.0129	2.4500	2.33
Cs L $\beta$ 4	2	5	2.0183	2.4565	4.65

S	SK $\alpha$ 4	1	45	2.0187	2.4569	2.32
	Au L $\beta$ 3	5	6	2.0206	2.4593	11.61
S	SK $\alpha$ 3	1	45	2.0212	2.4600	2.32
	Nb L $\beta$ 4	1	3	2.0230	2.4622	2.32
	Au L $\beta$ 2	5	20	2.0251	2.4647	11.58
S	SK $\alpha$ '	1	.3	2.0259	2.4657	2.32
	Pt M3-N4	1	1	2.0273	2.4675	2.31
	W M2-N4	1	.1	2.0273	2.4675	2.31
	Hg L $\beta$ 4	5	4	2.0289	2.4693	11.56
	Nb L $\beta$ 6	1	3	2.0289	2.4694	2.31
	Ga K $\alpha$ 1,2	4	150	2.0305	2.4714	9.24
	Cs L $\beta$ 1	2	50	2.0312	2.4722	4.62
	Co K $\alpha$ 1,2	3	150	2.0325	2.4738	6.92
S	K $\alpha$ 1,2	1	150	2.0334	2.4749	2.31
	Hg L $\delta$ 1	6	10	2.0355	2.4775	13.83
	Zr L $\delta$ 1	1	1	2.0376	2.4800	2.30
	Au L $\delta$ 3	6	2	2.0386	2.4813	13.81
	Sb L $\delta$ 2,3	2	2	2.0400	2.4829	4.60
	Pb L1	4	3	2.0434	2.4871	9.18
	Ir La1	4	100	2.0455	2.4896	9.17
	Mo La1	1	.00	2.0461	2.4903	2.29
	Mo La2	1	10	2.0490	2.4939	2.29
	Au L $\beta$ 1	5	50	2.0503	2.4954	11.44
	Te L $\delta$ 1	2	9	2.0530	2.4987	4.57
	Hg M $\beta$	1	50	2.0556	2.5019	2.28
	Ir La2	4	10	2.0625	2.5103	9.10
	Tl Ma1	1	100	2.0663	2.5149	2.27
	Tl Ma2	1	100	2.0708	2.5204	2.26
	Nb L $\beta$ 1	1	45	2.0785	2.5298	2.26
	W L $\delta$ 1	5	10	2.0787	2.5300	11.28

Zr L <sub>δ5</sub>	1	.1	2.0805	25323	2.26
Tl K <sub>α1,2</sub>	2	150	2.0812	2.5331	4.51
Ir M <sub>δ</sub>	1	1	2.0814	2.5333	2.25
Pt L <sub>β2</sub>	5	20	2.0852	2.5379	11.25
Pt L <sub>β3</sub>	5	6	2.0882	2.5416	11.23
Se K <sub>α1,2</sub>	5	150	2.0930	2.5474	11.21
Au L <sub>β4</sub>	5	4	2.0937	2.5483	11.20
Bi M3-N1	1	1	2.0954	2.5504	2.24
Tl L1	4	3	2.0962	2.5513	8.95
Ir M3-M4	1	.1	2.0966	2.5517	2.24
Ba L <sub>α1</sub>	2	100	2.1011	2.5572	4.47

NOTES: Sulfur appears to be a difficult element to analyze for depending on the mineral matrix. Of particular importance are Pb<sub>M<sub>α1</sub></sub>, Pb<sub>L<sub>α2</sub></sub>,


Ga<sub>K<sub>α1,2</sub></sub>, Cs<sub>L<sub>β1</sub></sub>, Co<sub>K<sub>α1,2</sub></sub>, Ir<sub>L1</sub>, Mo<sub>L<sub>α1</sub></sub> and Tl<sub>M<sub>α1</sub></sub>. Of these, Pb, Mo and Tl are first order X-ray lines of almost identical energy as S<sub>K<sub>α1,2</sub></sub> (2.31 KeV) and therefore not separable from sulfur with a PHA. Other interfering lines can be discriminated against with a PHA.



ORIGINAL X-RAY  
OF POOR QUALITY




70

FEET  
LEFT




S Co 99.9

ADP 300 cps.



S troilite  
PET 300 cps.  
S 36.47 wt. %



CHLORINE

TABLE 17

ELEMENT: C<sub>2</sub>

ANALYTICAL LINE: K<sub>α1,2</sub>

CRYSTAL: PET

SPECTROMETER SETTING: PET 2.1781

BACKGROUND SETTING: ±0.030

ELEMENTAL SCANS:


<u>MINERAL</u>	<u>WT. % C<sub>2</sub></u>	<u>BACKGROUND SETTING</u>
Scapolite	2.19	±0.030

## INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>ADP (101)</u>	<u>PET (002)</u>	<u>KeV</u>
	Ho L <sub>β10</sub>	3	.01	1.7582	2.1399	8.00
	Nb L <sub>δ2,3</sub>	1	.5	1.7613	2.1437	2.66
	Er L <sub>β3</sub>	3	6	1.7729	2.1578	7.94
	Gd L <sub>δ6</sub>	3	.01	1.7761	2.1617	7.92
C <sub>2</sub>	SK <sub>α4</sub>	1	4	1.7770	2.1628	2.64
	Sr K <sub>β1</sub>	6	16	1.7777	2.1637	15.83
	Ce L <sub>β4</sub>	2	5	1.7784	2.1646	5.28
	Sr K <sub>β3</sub>	6	8	1.7789	2.1652	15.82
	Ho L <sub>β2,15</sub>	3	20	1.7792	2.1654	7.91
	Er L <sub>β6</sub>	3	.1	1.7796	2.1660	7.91
	As K <sub>α1,2</sub>	4	150	1.7820	2.1689	10.53
C <sub>2</sub>	SK <sub>α3</sub>	1	4	1.7828	2.1699	2.63
	Gd L <sub>δ8</sub>	3	.1	1.7832	2.1704	7.89
	Ce L <sub>β1</sub>	2	50	1.7833	2.1705	5.26
C <sub>2</sub>	K <sub>α1,2</sub>	1	150	1.7896	2.1781	2.62
	Tb L <sub>δ5</sub>	3	.1	1.7923	2.1815	7.85
	Nd L <sub>α1</sub>	2	100	1.7941	2.1836	5.23
	Th L <sub>β4</sub>	6	4	1.7996	2.1904	15.64
	La L <sub>β6</sub>	2	.1	1.8006	2.1916	5.21
	Th L <sub>β2</sub>	6	20	1.8018	2.1930	15.62
	Nd L <sub>α2</sub>	2	10	1.8019	2.1931	5.21


Er L $\beta$ 1	3	50	1.8021	2.1933	7.81
Dy L $\beta$ 5	3	.1	1.8033	2.1948	7.80
Yb L2-M2	3	.01	1.8034	2.1949	7.80
Eu L $\delta$ 3	3	1	1.8055	2.1975	7.80
Gd L $\delta$ 1	3	5	1.8079	2.2004	7.78
Th La1	5	100	1.8089	2.2017	12.97

NOTES: C $\alpha$  is commonly found in sodalite, apatite and scapolite associated with REE, Sr, Th, As, V and F. Of these elements the REE could create most of the interference problems but they commonly occur in sufficiently low concentrations so as not to interfere with the C $\alpha$  determination. All but the Nb L $\delta$ 2,3 line could be discriminated against with a PHA window.




POTASSIUM











K kaersutite  
ADP 300 cps.  
K<sub>2</sub>O 1.59 wt. %





CALCIUM

TABLE 20

ELEMENT: Ca  
 ANALYTICAL LINE:  $K\alpha_{1,2}$   
 SPECTROMETER SETTING: 1.2714(ADP); 3.3595(LiF); 1.5474(PET)  
 BACKGROUND SETTING:  $\pm 0.046$ (ADP);  $\pm 0.052$ (LiF);  $\pm 0.034$ (PET)  
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT.% OXIDE</u>	<u>BACKGROUND SETTING</u>
CPX CP19	18.13	$\pm 0.046$ ADP $\pm 0.052$ LiF $\pm 0.034$ PET
Kaersutite	10.12	$\pm 0.034$ PET

## INTERFERENCES:

<u>Line</u>	<u>N</u>	<u>I</u>	<u>ADP (101)</u>	<u>LiF (200)</u>	<u>PET (002)</u>	<u>KeV</u>
Sm $L\delta 3$	2	1	1.2534	3.3120	1.5255	7.49
Eu $L\delta 1$	2	5	1.2545	3.3148	1.5263	7.48
Tb $L\beta 7$	2	.1	1.2553	3.3170	1.5278	7.47
Ni $K\alpha_{1,2}$	2	150	1.2558	3.3184	1.5285	7.47
Ho $L\beta 4$	2	5	1.2560	3.3190	1.5287	7.47
Sm $L\delta 2$	2	1	1.2568	3.3209	1.5296	7.47
Y $K\alpha_{1,2}$	4	150	1.2568	3.3210	1.5297	14.93
Nb $K\beta 3$	5	7	1.2608	3.3317	1.5346	18.6J
Tb $L\beta 10$	2	.01	1.2619	3.3346	1.5359	7.44
Ca $SK\alpha 4$	1	2	1.2628	3.3369	1.5370	3.72
Ca $SK\alpha 3$	1	2	1.2640	3.3401	1.5385	3.71
Po $L\alpha 1$	3	100	1.2646	3.3416	1.5391	11.13
Yb $L\alpha 1$	2	100	1.2654	3.3438	1.5402	7.41
Ca $K\alpha_{1,2}$	1	150	1.2714	3.3595	1.5474	3.69
Dy $L\beta 6$	2	.1	1.2732	3.3643	1.5496	7.37
Dy $L\beta 3$	2	6	1.2732	3.3644	1.5497	7.37
Yb $L\alpha 2$	2	10	1.2737	3.3657	1.5503	7.37
Tb $L\beta 2,15$	2	20	1.2738	3.3650	1.5504	7.37

Po La2	3	10	1.2778	3.3764	1.5552	11.01
Sm L68	2	.1	1.2915	3.4127	1.5719	7.27
Eu L65	2	.1	1.2931	3.4170	1.5739	7.26

NOTES: No interference problems typically exist due to the low concentrations of REE in most silicates. A high Ni concentration might present some problem but a PHA should be helpful in that case. If K is present and the analyzing crystal is PET, careful background selection is required. (See kaersutite scan.)




Ca Clinopyroxene CP19  
ADP 300 cps.  
CaO 18.13 wt. %

Ca<sub>2</sub>Al<sub>2</sub>O<sub>7</sub>


45


FEET  
LEFT






Ca Clinopyroxene CP19  
LIF 300 cps.  
CaO 18.13 wt. %





ORIGINAL PAGE IS  
OF POOR QUALITY



TITANIUM

TABLE 22

ELEMENT: Ti  
 ANALYTICAL LINE: K $\alpha$ 1,2  
 CRYSTAL: LiF, PET  
 SPECTROMETER SETTING: 2.7497, 1.2665  
 BACKGROUND SETTING:  $\pm 0.058$   
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT. % OXIDE</u>	<u>BACKGROUND SETTING</u>
rutile	100.00	$\pm 0.080$
ilmenite	50.04	$\pm 0.058$
kaersutite	5.64	$\pm 0.058$
Ti metal (99.9)		

## INTERFERENCES:

<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF</u> <u>(200)</u>	<u>PET</u> <u>(002)</u>	<u>KeV</u>
Cs L $\beta$ 4	1	5	2.6730	1.2312	9.28
La L $\alpha$ 2	1	10	2.6753	1.2323	4.63
Nd L1	1	2	2.6760	1.2326	4.63
Cs L $\beta$ 1	1	50	2.6837	1.2361	4.62
U L $\alpha$ 1	3	100	2.7319	1.2583	13.61
La Ln	1	1	2.7400	1.2621	4.52
Lu L $\beta$ 15	2	20	2.7430	1.2634	9.04
Hf L $\beta$ 1	2	50	2.7482	1.2658	9.02
Ti K $\alpha$ 1,2	1	150	2.7497	1.2665	4.51
Po L $\beta$ 1	3	50	2.7650	1.2740	13.44
Ba L $\alpha$ 1	1	100	2.7760	1.2786	4.47
Sc K $\beta$ 1,3	1	20	2.7796	1.2803	4.46
Rb K $\alpha$ 1,2	3	150	2.7808	1.2808	13.37
Pr L1	1	2	2.7841	1.2824	4.45
Cu K $\beta$ 1,3	2	20	2.7844	1.2825	8.90
Cu K $\beta$ 3	2	6	2.7852	1.2829	8.90
Ba L $\alpha$ 2	1	10	2.7855	1.2830	4.45
Zr K $\beta$ 1	4	18	2.8069	1.2929	17.67
U M $_2$ -N $_4$	1	5	2.8170	1.2975	4.40


NOTES: Commonly no interferences occur with a Ti determination with the exception of a high Ba matrix, i.e., benitoite. Other interferences will probably not be encountered. Ba interference cannot be separated with a PHA.



08

FEET  
LEFT


T1 T10<sub>2</sub>  
LIF 300 cps.  
T10<sub>2</sub> 100.0 WT. %



55

FEET  
LEFT

T1 ilmenite, Elliot Co.  
LiF 300 cps.  
TiO<sub>2</sub> 50.04 WT. %







Ti Kaersutite, Spring Mtn.

L1F 300 cps.

TiO<sub>2</sub> 5.64 wt. %



32

FEET  
LEFT

VANADIUM

TABLE 23


ELEMENT: V  
 ANALYTICAL LINE: K $\alpha$ 1,2  
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 2.5048  
 BACKGROUND SETTING:  $\pm 0.057$ ,  $\pm 0.029$   
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT. %</u>	<u>BACKGROUND SETTING</u>
V <sub>2</sub> O <sub>5</sub>	99.99	$\pm 0.057$
magnetite	?	$\pm 0.029$
ilmenite	?	-0.029
V, Ti overlay		

## INTERFERENCES:

<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
La L $\beta$ 4	1	5	2.4493	5.06
La L $\beta$ 1	1	50	2.4589	5.04
Pr L $\alpha$ 1	1	100	2.4630	5.03
Sm L1	1	2	2.4823	4.99
Ba L $\beta$ 6	1	.1	2.4826	4.99
Y K $\alpha$ 1,2	3	150	2.4907	14.93
Ti K $\beta$ 5	1	.02	2.4985	4.96
V K $\alpha$ 1,2	1	150	2.5048	4.95
Ge	2	150	2.5107	9.87
Pr LN	1	1	2.5120	4.94
Ti K $\beta$ 1,3	1	20	2.5139	4.93
Ba L $\beta$ 3	1	6	2.5164	4.93
Ba L $\beta$ 4	1	5	2.5553	4.85
Ce L $\alpha$ 1	1	100	2.5615	4.84


NOTES: The most likely problem to occur is the Ti K $\beta$ 1,3 interference on VK $\alpha$ 1,2. A PHA cannot be used as a discriminator due to the similar energies of the two lines. Using a V free sample, i.e. TiO<sub>2</sub>, ilmenite, etc. determine the intensity contribution at the V peak position and subtract off the appropriate amount from the V counts. See Geissman and Essene, "A graphic method of resolving X-ray interference on the electron microprobe", 13th Proceedings of the Microbeam Analysis Society. Ge interference can be discriminated against using a PHA.




ORIGINAL PAGE IS  
OF POOR QUALITY





V magnetite, section 3A  
LiF 300 cps.  
 $V_2O_3$  ? wt. %











CHROMIUM

TABLE 24

ELEMENT: Cr

ANALYTICAL LINE: K $\alpha$ 1,2

CRYSTAL: LiF

SPECTROMETER SETTING: 2.291

BACKGROUND SETTING:  $\pm 0.064$

## ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT % OXIDE</u>	<u>BACKGROUND SETTING</u>
Chromite	38.76	$\pm 0.064$
Armcolite	6.69	$\pm 0.046$


## INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
	Ba L $\gamma$ 1	1	5	2.2415	5.53
	Pr L $\beta$ 4	1	5	2.2550	5.50
	Pr L $\beta$ 1	1	50	2.2588	5.49
	V K $\beta$ 5	1	.02	2.2695	5.49
	V K $\beta$ 1,3	1	20	2.2844	5.43
Cr	K $\alpha$ 1,2	1	150	2.2910	5.41
	La L $\beta$ 2,15	1	20	2.3030	5.38
	Ba L $\gamma$ 5	1	.1	2.3085	5.37
	Ce L $\beta$ 3	1	6	2.3109	5.36
	Ce L $\beta$ 4	1	5	2.3497	5.28
	Ce L $\beta$ 1	1	50	2.3561	5.26

NOTES: Chromium in a high concentration of V, Ba, or REE may be a difficult determination because of the similar energies of all the interfering X-ray lines. The use of a PHA as a discriminator is hopeless.


FEET  
LEFT

Cr chromite  
LiF 300 cps.  
Cr<sub>2</sub>O<sub>3</sub> 38.76 wt. %





58

FEET  
LEFT

MANGANESE



TABLE 25

ELEMENT: Mn  
 ANALYTICAL LINE:  $K\alpha_{1,2}$   
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 2.1031  
 BACKGROUND SETTING:  $\pm 0.080$  for Rhodonite;  $\pm 0.030$  low concentration samples  
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT. % OXIDE</u>	<u>BACKGROUND SETTINGS</u>	
Rhodonite	40.77	$\pm 0.080$	
Diopside	1.41	$\pm 0.035$	
Chromite	0.17	Determination	difficult

## INTERFERENCES


	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
	Nd $L\beta_{2,15}$	1	20	2.0360	6.09
	Gd $L\alpha_1$	1	100	2.0468	6.06
	Ce $L\gamma_1$	1	5	2.0487	6.05
	Gd $L\alpha_2$	1	10	2.0578	6.02
	Cr $K\beta_{1,3}$	1	18	2.0849	5.95
Mn	$K\alpha_{1,2}$	1	150	2.1031	5.89
	Eu $L\alpha_1$	1	100	2.1209	5.85
	Nd $L\beta_3$	1	6	2.1268	5.83
	Eu $L\alpha_2$	1	10	2.1315	5.82
	Mo $K\alpha_{1,2}$	3	150	2.1322	17.44
	La $L\gamma_1$	1	5	2.1418	5.79
	Nd $L\beta_4$	1	5	2.1669	5.72
	Nd $L\beta_1$	1	50	2.1669	5.72

NOTES: Typically no interference problems except with a Cr-rich matrix, i.e., chromite, where  $CrK\beta_{1,3}$  will interfere. PHA window is of no help.

53

FEET  
LEFT

Mn rhodonite  
LIF 300 cps.  
MnO 40.77 wt. %





ORIGINAL FILED  
IN 100-100-100

Mn diopside glass doped with Cr, Mn, Ni.

LiF 300 cps.

MnO 1.41 wt. %





IRON

c-2

TABLE 26

ELEMENT: Fe  
 ANALYTICAL LINE: K $\alpha$ 1,2  
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 1.937  
 BACKGROUND SETTING:  $\pm 0.062$   
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT.% OXIDE</u>	<u>BACKGROUND SETTING</u>
CP19	8.49	$\pm 0.062$
Rhodonite	3.60	$\pm 0.062$

## INTERFERENCES:

	<u>Line</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
Nd	L $\gamma$ 1	1	5	1.8779	6.60
Pt	LV	2	.01	1.8786	13.20
Pt	L $\gamma$ 8	2	.1	1.8822	13.17
Pt	L $\gamma$ 1-N $\gamma$ 1	2	.01	1.8845	13.16
Eu	L $\beta$ 3	1	6	1.8867	6.57
Mo	K $\beta$ 1	3	17	1.8969	19.61
Mn	K $\beta$ 5	1	.03	1.8971	6.53
Mo	K $\beta$ 3	3	7	1.8986	19.59
DY	L $\alpha$ 1	1	100	1.9088	6.49
Mn	K $\beta$ 1,3	1	20	1.9102	6.49
Pt	L $\gamma$ 1	2	10	1.9159	12.94
Eu	L $\beta$ 1	1	50	1.9203	6.46
Eu	L $\beta$ 4	1	5	1.9255	6.44
Fe	K $\alpha$ 1,2	1	150	1.9374	6.40
Pt	L2-N3	2	.01	1.9443	12.75
Pt	L2-N2	2	.01	1.9584	12.66
Pr	L $\gamma$ 1	1	5	1.9611	6.32

Sm	L $\beta$ 3	1	6	1.9624	6.32
Pt	L $\gamma$ 5	2	.1	1.9754	12.55
Tb	L $\alpha$ 2	1	10	1.9875	6.24
Sm	L $\beta$ 1	1	50	1.9981	6.20
Sm	L $\beta$ 4	1	5	2.0010	6.20

NOTES: Typically no problem except in a high Mn matrix. See rhodonite scan.  
Platinum lines are included for the benefit of experimental petrologists  
using Pt capsules.




61


FEET  
LEFT

Fe clinopyroxene CP19

LiF 300 cps.

FeO 8.49 wt. %





ORIGINAL PAGE 1  
OF 2 (10/1/77)

COBALT

TABLE 27

ELEMENT: Co  
 ANALYTICAL LINE:  $K\alpha_{1,2}$   
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 1.7903  
 BACKGROUND SETTING:  $\pm 0.080$   
 ELEMENTAL SCANS:


<u>MINERAL</u>	<u>WT. %</u>	<u>BACKGROUND SETTING</u>
Kamacite	—	$\pm 0.080$
Co	99.9	$\pm 0.080$

## INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
Dy	$L\beta_1$	1	50	1.7106	7.25
Sm	$L\gamma_1$	1	5	1.7272	7.18
Fe	$K\beta_5$	1	.03	1.7442	7.11
Sr	$K\alpha_{1,2}$	2	150	1.7533	14.14
Fe	$K\beta_{1,3}$	1	20	1.7566	7.06
Tb	$L\beta_1$	1	50	1.7768	6.98
Er	$L\alpha_1$	1	100	1.7843	6.95
Tb	$L\beta_4$	1	5	1.7864	6.94
Co	$K\alpha_{1,2}$	1	150	1.7903	6.92
Eu	$L\beta_{2,15}$	1	20	1.8118	6.84


NOTES: Typically no problem except in a high Fe matrix. (See kamacite scan.)  
 Undoubtedly there is some contribution to Co  $K\alpha_{1,2}$  from Fe  $K\beta_{1,3}$ .  
 Due to similar energies of the X-ray lines a PHA is not helpful.  
 Best method of analysis would be to choose a standard with as similar composition as the sample as possible.

Co Kamacite  
Yamato Meteorite  
LIF 300 cps.





Co 99.9  
LIF 300 cps.



NICKEL



TABLE 28

ELEMENT: Ni  
 ANALYTICAL LINE: K $\alpha$ 1,2  
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 1.6592  
 BACKGROUND SETTING:  $\pm 0.0065$   
 ELEMENTAL SCANS:


MINERAL	WT. %	BACKGROUND SETTING
Ni (99.9)	—	$\pm 0.065$
Fe85 Ni10 Co5	10.0 Ni	—
FeNi phase (Yamato Meteorite)	—	$\pm 0.029$


## INTERFERENCES:

	LINE	N	I	LiF (200)	KeV
Er	L $\beta$ 4	1	5	1.6007	7.74
Co	K $\beta$ 5	1	.03	1.6089	7.71
Lu	L $\alpha$ 1	1	100	1.6195	7.65
Mo	L $\beta$ 3	1	6	1.6203	7.65
Co	K $\beta$ 1,3	1	20	1.6208	7.65
DY	L $\beta$ 2,15	1	20	1.6237	7.63
Lu	L $\alpha$ 2	1	10	1.6303	7.60
Eu	L $\gamma$ 1	1	5	1.6574	7.48
Ni	K $\alpha$ 1,2	1	150	1.6592	7.47
Ho	L $\beta$ 4	1	5	1.6595	7.47
Y	K $\alpha$ 1,2	2	150	1.5605	14.93
Yb	L $\alpha$ 1	1	100	1.6719	7.41

DY	Lβ3	1	6	1.6822	7.37
Ag	Kα1,2	3	150	1.6826	22.10
Yb	Lα2	1	10	1.6829	7.37
DY	Lβ1	1	50	1.7106	7.25
DY	Lβ4	1	5	1.7210	7.20
Tm	Lα1	1	100	1.7268	7.18


NOTES: Typically no problem although a high Co matrix may present an overlap problem due to tailing effects, but at the 5% Co level there does not appear to be a problem.





ORIGINAL PAGE IS  
OF POOR QUALITY

Ni Yamato Meteorite  
LiF 300 cps.



STRONTIUM



TABLE 38

ELEMENT: Sr  
 ANALYTICAL LINE: La1  
 CRYSTAL: ADP; PET  
 SPECTROMETER SETTING: 2.5972; 3.1610  
 BACKGROUND SETTING:  
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT. % OXIDE</u>	<u>BACKGROUND SETTING</u>
Lake County Plagioclase	600 ppm Sr	
Celestite	56.41	±0.046
Glass X (Probe Society)	0.80	


## INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>ADP (101)</u>	<u>PET (002)</u>	<u>KeV</u>
	Ca Ka1,2	2	150	2.5427	3.0948	3.69
	Re Ma	1	100	2.5465	3.0994	1.84
	Si Kb	1	2	2.5556	3.1105	1.84
	W Mb	1	45	2.5571	3.1123	1.84
	Lu My	1	1	2.5613	3.1174	1.83
	Si Kb1	1	2	2.5651	3.1220	1.83
	Rb Lb3	1	3	2.5687	3.1264	1.83
	Pb MZ2	1	.1	2.5741	3.1330	1.82
	Si SkB'	1	.1	2.5793	3.1394	1.82
	Rb Lb4	1	3	2.5812	3.1416	1.82
	V Kb1,3	3	20	2.5935	3.1566	5.43
Sr	La1	1	100	2.5972	3.1610	1.81
	Cr Ka1,2	3	150	2.6010	3.1657	5.41
	Sb La1	2	100	2.6032	3.1684	3.60
	K Kb1,3	2	15	2.6142	3.1818	3.59
	Zr LI	1	3	2.6183	3.1867	1.79
	W Ma1	1	100	2.6426	3.2164	1.78

Rb L $\beta$ 6	1	3	2.6431	3.2169	1.78
W M $\alpha$ 2	1	100	2.6461	3.2205	1.77
Si SK $\alpha$ 6	1	.1	2.6566	3.2334	1.77
Ta M $\beta$	1	45	2.6578	3.2348	1.77
Yb M $\gamma$	1	1	2.6582	3.2353	1.77
Fe K $\beta$ 1,3	4	20	2.6591	3.2364	7.06
Si SK $\alpha$ 5	1	.2	2.6603	3.2378	1.76
Si SK $\alpha$ 4	1	3	2.6743	3.2549	1.75
Rb L $\beta$ 1	1	45	2.6778	3.2592	1.75
Si SK $\alpha$ 3	1	6	2.6781	3.2596	1.75
Si SK $\alpha$ 1	1	.5	2.6848	3.2677	1.75
Si SK $\alpha$ 1,2	1	150	2.6969	3.2824	1.74

NOTES: Sr occurs at the trace level in silicates such as feldspars and clinopyroxenes and rarely occurs as a major constituent except in some SO<sub>4</sub> and CO<sub>3</sub> phases. At trace levels the adjacent Si SK line causes a background as well as an overlap problem and background position should be chosen on the high angle side only. The energy of the Si SK line is identical to the Sr L $\alpha$  line and therefore use of a PHA is hopeless. Better PK/BKg ratio can be achieved by using an electron beam of 30KV rather than 15KV acceleration potential.


57

FEET  
LEFT

58

FEET  
LEFT


Sr CELESTITE  
ADP 300 cps.  
SrO 56.41 wt. %



Sr GLASS X (Probe Society)

ADP 100 cps.

SrO 0.80



ORIGINAL PAGE IS  
OF POOR QUALITY

ZIRCONIUM



TABLE 40

ELEMENT: Zr  
 ANALYTICAL LINE: Zr La1,2  
 CRYSTAL: PET, ADP  
 SPECTROMETER SETTING: 2.7961; 2.2973  
 BACKGROUND SETTING:  $\pm 0.060$   
 ELEMENTAL SCANS:


<u>MINERAL</u>	<u>WT.% OXIDE</u>	<u>BACKGROUND SETTING</u>
zircon	67.27	$\pm 0.060$
armalcolite	6.31	$\pm 0.030$

## INTERFERENCES:


	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>PET (002)</u>	<u>ADP (101)</u>	<u>KeV</u>
	Y L $\beta$ 3	1	3	2.7559	2.2643	2.07
	Pt M $\gamma$ -O $\gamma$ 3	1	.01	2.7576	2.2657	2.07
	Sm L $\beta$ 1	3	50	2.2684	2.7609	6.20
	Hf K $\beta$ 1,3	4	20	2.7639	2.2708	8.26
	Y L $\beta$ 4	1	3	2.7722	2.2777	2.06
	Rb L $\gamma$ 2,3	1	.1	2.7847	2.2880	2.05
	Pt M $\alpha$ 1	1	100	2.7853	2.2884	2.05
	Er L $\beta$ 2,15	4	20	2.7894	2.2918	8.19
	Pt M $\alpha$ 2	1	100	2.7903	2.2926	2.05
	P SK $\alpha$ 6	1	.01	2.7927	2.2946	2.04
	Sc K $\alpha$ 1,2	2	150	2.7931	2.2949	4.09
Zr	La1	1	100	2.7961	2.2973	2.04
	P SK $\alpha$ 5	1	.01	2.7983	2.2991	2.04
Zr	La2	1	10	2.7995	2.3001	2.04
	Ta La1	4	100	2.8041	2.3039	8.15
	Y L $\beta$ 6	1	3	2.8070	2.3063	2.03
	Nd L $\beta$ 2,15	3	20	2.8134	2.3115	6.09

P SK $\alpha$ 4	1	5	2.8137	2.3118	2.03
P SK $\alpha$ 3	1	5	2.8173	2.3148	2.03
Tm L $\beta$ 1	4	50	2.8196	2.3167	8.10
P SK $\alpha$ '	1	.5	2.8239	2.3202	2.02
Gd L $\alpha$ 1	3	100	2.8283	2.3238	6.06
P K $\alpha$ 1,2	1	150	2.8364	2.3304	2.01
Ca K $\beta$ 1,3	2	15	2.8463	2.3385	4.01

NOTES: Chief source of interference with Zr determination is P SK $\alpha$ 5. Although of low intensity it occurs at the same position and is the same energy as the Zr L $\alpha$ 1 therefore a PHA is not helpfull. Aside from P the REE Sm, Er, Nd, Tm and Gd all have x-ray peaks within the background region but can be discriminated against with a PHA. Ca K $\beta$ 1,3 is probably far enough away so as not to present a problem.



Zr Armalcolite  
ADP 300 cps.  
ZrO<sub>2</sub> 6.31 wt. %



TEXAS INSTRUMENTS INCORPORATED HOUSTON, TEXAS U.S.

ORIGINAL PAGE IS  
OF POOR QUALITY

BAR1011

TABLE 56

ELEMENT: Ba

ANALYTICAL LINE:  $\text{La1}$ CRYSTAL:  $\text{LiF (200)}$ 

SPECTROMETER SETTING: 2.7760

BACKGROUND SETTING:  $\pm 0.056$ 

ELEMENTAL SCANS:


<u>MINERAL</u>	<u>WT % OXIDE</u>	<u>BACKGROUND SETTING</u>
Barite	65.69	$\pm 0.056$
Benitoite		

INTERFERENCE:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
Er	$\text{La1}$	2	5	2.7282	9.09
Ti	$\text{Ka1,2}$	1	150	2.7497	4.51
Zr	$\text{KB2}$	4	4	2.7597	17.97
Ba	$\text{La1}$	1	100	2.7760	4.47
Rb	$\text{Ka1,2}$	3	150	2.7808	13.37
Cu	$\text{KB1,3}$	2	20	2.7844	8.90
Cu	$\text{KB3}$	2	6	2.7852	8.90
Ba	$\text{La2}$	1	10	2.7855	4.45
Zr	$\text{KB1}$	4	18	2.8069	17.67
Zr	$\text{KB3}$	4	9	2.8091	17.65

NOTES: Barium typically occurs as a minor to trace level constituent in feldspars and as a major component of some carbonates, sulfates and halides. There should be little or no interference on the  $\text{BaLa1}$  line with the exception of  $\text{TiKa1,2}$  which is of such similar energy (see TABLE) that use of a PHA will not discriminate against it.






ORIGINAL PAGE IS  
OF POOR QUALITY

70

LEFT



LANTHANUM

TABLE 57

ELEMENT: La  
 ANALYTICAL LINE: La1  
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 2.6657  
 BACKGROUND SETTING:  
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT % OXIDE</u>	<u>BACKGROUND SETTING</u>
Glass REE 3	4.28	±0.027
Perovskite	?	

## INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
	Tm L <sub>γ</sub> 1	2	5	2.6306	9.42
	Ho L <sub>γ</sub> 1-04,5	2	.01	2.6416	9.39
	Er L <sub>γ</sub> 2	2	1	2.6420	9.38
	Ho L <sub>γ</sub> 4	2	.1	2.6450	9.37
	Pt L <sub>α</sub> 2	2	10	2.6486	9.36
	Hf L <sub>β</sub> 2	2	20	2.6528	9.35
	Hf L <sub>β</sub> 15	2	1	2.6557	9.34
La	La1	1	100	2.6657	4.65
	Cs L <sub>β</sub> 4	1	5	2.6666	4.65
	Lu L <sub>β</sub> 9	2	.01	2.6716	9.28
La	La2	1	10	2.6753	4.63
	Nd L1	1	2	2.6760	4.63
	Er L <sub>γ</sub> 6	2	.01	2.6794	9.25
	Lu L <sub>β</sub> 5	2	.1	2.6837	9.24
	Cs L <sub>β</sub> 1	1	50	2.6837	4.62
	Lu L <sub>β</sub> 10	2	.01	2.6860	9.23
	Lu L <sub>3</sub> -02,3	2	.01	2.6905	9.22

NOTES: Chief interferences come from 1st order  $CsL\beta_4$  and  $NdL_1$  x-ray lines. Due to the similar energies of these lines with  $La_{L\alpha}$ , a PHA is of no use as an energy discriminator. All 2nd order lines can be dealt with using a PHA.

CERIUM



TABLE 58

ELEMENT: Ce  
 ANALYTICAL LINE:  $L_{\alpha 1}$   
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 2.5615  
 BACKGROUND SETTING:  
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT % OXIDE</u>	<u>BACKGROUND SETTING</u>
Glass REE 3	4.00	$\pm 0.027$
perovskite	?	

## INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
	Tm $L_{\gamma 3}$	2	1	2.5356	9.78
	Tm $L_{\gamma 2}$	2	1	2.5484	9.73
	Er $L_{\gamma 4}$	2	.1	2.5504	9.72
	Ba $L_{\beta 4}$	1	5	2.5553	4.85
Ce	$L_{\alpha 1}$	1	100	2.5615	4.84
	Po $L_I$	2	3	2.5658	9.66
	Ba $L_{\beta 1}$	1	50	2.5682	4.83
Ce	$L_{\alpha 2}$	1	10	2.5706	4.82
	Hf $L_{\beta 9}$	2	.01	2.5805	9.61
	Tm $L_{\gamma 6}$	2	.01	2.5810	9.61
	Zn $K_{\beta 1,3}$	2	20	2.5905	9.57

NOTES: Chief interference is from Ba  $L_{\beta 1}$  and  $L_{\beta 4}$ . Similar X-ray energies prevent energy discrimination with a PHA. All other lines can be discriminated against using the PHA.

PRASEODYMIUM

TABLE 59

ELEMENT: Pr  
 ANALYTICAL LINE:  $L_{\alpha 1}$   
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 2.4630  
 BACKGROUND SETTING:  
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT % OXIDE</u>	<u>BACKGROUND SETTING</u>
Glass REE 3	4.45	$\pm 0.027$

## INTERFERENCES:

<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
Hf $L_{\gamma 5}$	2	.1	2.4307	10.20
Th $M_1-O_3$	1	.01	2.4420	5.08
U $M_2-O_4$	1	.01	2.4430	5.07
Lu $L_{\gamma 1}$	2	5	2.4446	10.14
Yb $L_{\gamma 3}$	2	1	2.4446	10.14
La $L_{\beta 4}$	1	5	2.4493	5.06
Tm $L_1 O_{4,5}$	2	.01	2.4526	10.11
Yb $L_{\gamma 2}$	2	1	2.4576	10.09
Tm $L_{\gamma 4}$	2	.1	2.4588	10.08
La $L_{\beta 1}$	1	50	2.4589	5.04
Pr $L_{\alpha 1}$	1	100	2.4630	5.03
Pr $L_{\alpha 2}$	1	10	2.4729	5.01
Sm $L_1$	1	2	2.4823	4.99

Ba L <sub>β</sub> 6	1	.1	2.4826	4.99
Yb L <sub>γ</sub> 6	2	.01	2.4854	9.98
Pt Ln	2	1	2.4858	9.97
Yb L <sub>2</sub> -0 <sub>2,3</sub>	2	.01	2.4906	9.95

NOTES: Interferences which occur and which the analyst is helpless against are: ThM<sub>1</sub>-0<sub>3</sub>, U<sub>M</sub><sup>2</sup>-0<sub>4</sub>, La L<sub>β</sub>4, LaL<sub>β</sub>1, Sm L<sub>1</sub>, and BaL<sub>β</sub>6. La would be the most likely problem. All other lines can be discriminated against with a PHA.

NEODYMIUM

TABLE 60

ELEMENT: Nd  
 ANALYTICAL LINE: La1  
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 2.3704  
 BACKGROUND SETTING:  
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT % OXIDE</u>	<u>BACKGROUND SETTING</u>
Glass REE 2	4.26	±0.030
Perovskite		

## INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
	Cs L <sub>γ</sub> 1	1	5	2.3480	5.28
	Sr K <sub>β</sub> 1	3	16	2.3488	15.83
	Ce L <sub>β</sub> 4	1	5	2.3497	5.28
	Pb L <sub>α</sub> 1	2	100	2.3500	10.55
	Ce L <sub>β</sub> 1	1	50	2.3561	5.26
	Zr K <sub>α</sub> 1,2	3	150	2.3620	15.74
	Po L <sub>γ</sub> 1	3	10	2.3624	15.74
Nd	La1	1	100	2.3704	5.23
	La L <sub>β</sub> 6	1	.1	2.3790	5.21
	Ba L <sub>β</sub> 7	1	.1	2.3806	15.62
Nd	La2	1	10	2.3807	5.21
	Ba L <sub>β</sub> 10	1	.01	2.3869	5.18
	Eu L1	1	2	2.3948	5.18
	Ba L <sub>β</sub> 2,15	1	20	2.4044	5.16

NOTES: Intense X-ray lines from elements such as Zr, Pb are more energetic than Nd La and can be discriminated against with a PHA. Elements such as Ce, La and Ba, if present in the matrix, cannot be discriminated against and interfere with the measurement.

SALVARIUM



TABLE 62

ELEMENT: Sm  
 ANALYTICAL LINE:  $L\alpha_1$   
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 2.1998  
 BACKGROUND SETTING:  
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT % OXIDE</u>	<u>BACKGROUND SETTING</u>
Glass REE 2	4.26	-0.017

## INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
	Nd $L\beta_1$	1	50	2.1669	5.72
	Cs $L\gamma_4$	1	.1	2.1741	5.70
	Ce $L\beta_7$	1	.1	2.1806	5.69
	Ce $L\beta_9$	1	.01	2.1885	5.66
	Pt $L\beta_6$	1	.1	2.1906	5.66
	Ce $L\beta_{10}$	1	.01	2.1958	5.65
Sm	$L\alpha_1$	1	100	2.1998	5.64
	Pt $L\beta_2$	2	20	2.2040	11.25
	La $L\gamma_5$	1	.1	2.2056	5.62
	Ce $L\beta_{2,15}$	1	20	2.2087	5.61
Sm	$L\alpha_2$	1	10	2.2106	5.61
	Pr $L\beta_3$	1	6	2.2172	5.59
Sm	$L\alpha$	1	1	2.2182	5.59
	Po $L\alpha_1$	2	100	2.2277	11.13

NOTES: Ce, Cs, Pr, or La in the sample matrix can present a problem. Ce is the chief problem, but the continuum is loaded with other lines of similar energy. PHA is only useful for Pt and Po.

EUROPIUM

TABLE 63

ELEMENT: Eu  
 ANALYTICAL LINE: La1  
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 2.1209  
 BACKGROUND SETTING:  
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT % OXIDE</u>	<u>BACKGROUND SETTING</u>
Glass REE 1	4.21	±0.030

## INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
	Pr L $\beta$ 7	1	.1	2.0919	5.93
	Pr L $\beta$ 9	1	.1	2.1000	5.90
	Mn K $\alpha$ 1,2	1	150	2.1031	5.89
	Nd L $\beta$ 6	1	.1	2.1039	5.89
	Pr L $\beta$ 10	1	.01	2.1071	5.88
	Ce L $\gamma$ 5	1	.1	2.1103	5.87
	Pr L $\beta$ 2.15	1	20	2.1194	5.85
Eu	La1	1	100	2.1209	5.85
	Nd L $\beta$ 3	1	6	2.1268	5.83
Eu	LN	1	1	2.1315	5.82
Eu	La2	1	10	2.1315	5.82
	Ba L $\gamma$ 3	1	1	2.1342	5.81
	Ba L $\gamma$ 2	1	1	2.1387	5.80
	La L $\gamma$ 1	1	5	2.1418	5.79

NOTES: All of the x-ray lines occurring near Eu are of similar energies and therefore impossible to discriminate against with a PHA. Aside from La, Pr, Nd and Ba occurring very near the peak, Mn $K_{\alpha 1,2}$  is apt to be more of a problem, and depending on the Mn concentration level, the Eu might be sitting in the Mn peak tail.

GADOLINIUM

TABLE 64

ELEMENT: Gd  
 ANALYTICAL LINE: La1  
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 2.0468  
 BACKGROUND SETTING:  
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT % OXIDE</u>	<u>BACKGROUND SETTING</u>
Glass REE 1	4.47	±0.030

## INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
	Nd L $\beta$ 9	1	.01	2.0165	6.15
	Pr L $\gamma$ 5	1	.1	2.0205	6.14
	Nd L $\beta$ 10	1	.01	2.0237	6.13
	Ce L $\gamma$ 8	1	.1	2.0237	6.13
	Nd L $\beta$ 2,15	1	20	2.0360	6.09
	La L $\gamma$ 3	1	1	2.0410	6.07
	Pm L $\beta$ 3	1	6	2.0410	6.07
	La L $\gamma$ 2	1	1	2.0460	6.06
Gd	La1	1	100	2.0468	6.06
	Ce L $\gamma$ 1	1	5	2.0487	6.05
Gd	LN	1	1	2.0494	6.05
Gd	La2	1	10	2.0578	6.02
	Cr K $\beta$ 5	1	.03	2.0709	5.99
	Ba L $\gamma$ 4	1	.1	2.0756	5.97

NOTES: Chief interferences are Nd, Pm and Ce. A PHA will not help since all x-ray lines are of similar energy.

TERBIUM

TABLE 65

ELEMENT: Tb  
 ANALYTICAL LINE: La1  
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 1.9765  
 BACKGROUND SETTING:  
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT % OXIDE</u>	<u>BACKGROUND SETTING</u>
Glass REE 1	4.27	±0.026

## INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
	Sm L $\beta$ 6	1	.1	1.9464	6.37
	Ce L $\gamma$ 3	1	1	1.9553	6.34
	Pm L $\beta$ 2,15	1	20	1.9559	6.34
	Ce L $\gamma$ 2	1	1	1.9602	6.32
	Pr L $\gamma$ 1	1	5	1.9611	6.32
	Sm L $\beta$ 3	1	6	1.9624	6.32
Tb	LN	1	1	1.9730	6.28
Tb	La1	1	100	1.9765	6.27
	La L $\gamma$ 4	1	.1	1.9830	6.25
Tb	La2	1	10	1.9875	6.24
	Sm L $\beta$ 1	1	50	1.9981	6.20
	Sm L $\beta$ 4	1	5	2.0010	6.20

NOTES: Chief interferences for Tb measurement are Sm, Pr, and Pm. As with other REE, interfering x-ray lines are of sufficiently similar energy that a PHA is of no help.



DYSPROSIUM

TABLE 66

ELEMENT: Dy  
 ANALYTICAL LINE:  $L_{\alpha 1}$   
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 1.9088  
 BACKGROUND SETTING:  
 ELEMENTAL SCANS:

MINERAL  
 Glass REE 4

WT % OXIDE  
 4.36

BACKGROUND SETTING  
 $\pm 0.030$

## INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
	Nd $L_{\gamma 1}$	1	5	1.8779	6.60
	Pr $L_{\gamma 2}$	1	1	1.8791	6.60
	Sm $L_{\beta 2,15}$	1	20	1.8822	6.59
	Eu $L_{\beta 3}$	1	6	1.8867	6.57
	Yb $L_1$	1	2	1.8942	6.54
Dy	LN	1	1	1.8974	6.53
	Ce $L_{\gamma 4}$	1	.1	1.8991	6.53
Dy	$L_{\alpha 1}$	1	100	1.9088	6.49
	MnK $\beta 1,3$	1	20	1.9102	6.49
	Th $L_{\alpha 1}$	2	100	1.9120	12.97
Dy	$L_{\alpha 2}$	1	10	1.9199	6.46
	Eu $L_{\beta 1}$	1	50	1.9203	6.46
	Eu $L_{\beta 4}$	1	5	1.9255	6.44

Nd L $\gamma$ 5	1	.1	1.9355	6.41
Th La2	2	10	1.9358	12.81
Pr L $\gamma$ 8	1	.1	1.9362	6.40
Fe K $\alpha$ 1,2	1	150	1.9374	6.40

NOTES: Dy in a mafic mineral will be almost impossible to analyze due to the large Fe and Mn concentrations. In fact, the Dy (peak) is apt to be sitting in the Fe tail and MnK $\beta$ 1 is right on top of the Dy peak. All X-ray lines, with the exception of Th, are of similar energy and impossible to discriminate against with a PHA. These problems aside, Eu and Sm will also cause interference. Very difficult to analyze.

HOLMIUM

TABLE 67

ELEMENT: Ho  
 ANALYTICAL LINE: La1  
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 1.8450  
 BACKGROUND SETTING:  
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT % OXIDE</u>	<u>BACKGROUND SETTING</u>
Glass REE 4	4.41	±0.030

## INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
	Eu L $\beta$ 2,15	1	20	1.8118	6.84
	Gd L $\beta$ 3	1	6	1.8150	6.83
	Pr L $\gamma$ 4	1	.1	1.8193	6.81
	U La1	2	100	1.8213	13.61
Ho	LN	1	1	1.8264	6.79
	Yb Lt	1	.01	1.8309	6.77
	Lu L1	1	2	1.8360	6.75
Ho	La1	1	100	1.8450	6.72
	U La2	2	10	1.8451	13.44
	Gd L $\beta$ 1	1	50	1.8468	6.71
	Sm L $\beta$ 5	1	.1	1.8470	6.71
	Rb Ka1,2	2	150	1.8539	13.37
	Gd L $\beta$ 4	1	5	1.8540	6.69
	Nd L $\gamma$ 8	1	.1	1.8552	6.68

Ho	La2	1	10	1.8561	6.68
	Sm L $\beta$ 7	1	.1	1.8563	6.68
	Sm L $\beta$ 9	1	.01	1.8617	6.66
	Sm L $\beta$ 10	1	.01	1.8699	6.63
	Eu L $\beta$ 6	1	.1	1.8737	6.62
	Pr L $\gamma$ 3	1	1	1.8740	6.62

NOTES: Chief Ho interferences are Gd, Rb, Lu, U, and Eu. U and Rb are the only X-ray lines which can be discriminated against using a PHA.

ERBIUM



TABLE 68

ELEMENT: Er  
 ANALYTICAL LINE:  $L_{\alpha 1}$   
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 1.7843  
 BACKGROUND SETTING:  
 ELEMENTAL SCANS:

MINERAL  
 Glass REE 4

WT % OXIDE  
 4.36

BACKGROUND SETTING  
 $\pm 0.030$

## INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
	Sr $K_{\alpha 1,2}$	2	150	1.7533	14.14
	Fe $K_{\beta 1,3}$	1	20	1.7566	7.06
	Lu Lt	1	.01	1.7760	6.98
	Tb $L_{\beta 1}$	1	50	1.7768	6.98
	Eu $L_{\beta 5}$	1	.1	1.7772	6.98
	Sm $L_{\gamma 5}$	1	.1	1.7793	6.97
Er	$L_{\alpha 1}$	1	100	1.7843	6.95
	Eu $L_{\beta 7}$	1	.1	1.7851	6.94
	Tb $L_{\beta 4}$	1	5	1.7864	6.94
	Co $K_{\alpha 1,2}$	1	150	1.7903	6.92
	Eu $L_{\beta 9}$	1	.01	1.7916	6.92
Er	$L_{\alpha 2}$	1	10 <sup>-</sup>	1.7955	6.90

Nd L $\gamma$ 3	1	1	1.7964	6.90
Pm L $\gamma$ 1	1	5	1.7989	6.89
Eu L $\beta$ 10	1	.01	1.7993	6.89
Nd L $\gamma$ 2	1	1	1.8013	6.88
Gd L $\beta$ 6	1	.1	1.8054	6.87
Eu L $\beta$ 2,15	1	20	1.8118	6.84

NOTES: Major sources of interference are Co<sub>K $\alpha$ 1,2</sub>; Fe<sub>K $\beta$ 1,3</sub>; Tb<sub>L $\beta$ 4</sub> and Tb<sub>L $\beta$ 1</sub>. All lines except Sr are of similar energy and therefore a PHA cannot discriminate against them.

THULIUM

TABLE 69

ELEMENT: Tm  
 ANALYTICAL LINE: L $\alpha$ 1  
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 1.7268  
 BACKGROUND SETTING:  
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT % OXIDE</u>	<u>BACKGROUND SETTING</u>
Glass REE 1	4.36	-0.043

## INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
Tm	Ln	1	1	1.6963	7.31
	Sm L $\gamma$ 6	1	.01	1.6967	7.31
	Sm L $\gamma$ 8	1	.1	1.7063	7.27
	Eu L $\gamma$ 5	1	.1	1.7085	7.26
	Dy L $\beta$ 1	1	50	1.7106	7.25
	Gd L $\beta$ 5	1	.1	1.7130	7.24
	Gd L $\beta$ 7	1	.1	1.7203	7.21
	Dy L $\beta$ 4	i	5	1.7210	7.20
	Gd L $\beta$ 9	1	.01	1.7240	7.19
Tm	L $\alpha$ 1	1	100	1.7268	7.18
	Sm L $\gamma$ 1	1	5	1.7272	7.18
	Gd L $\beta$ 10	1	.01	1.7315	7.16
Tm	L $\alpha$ 2	1	10	1.7381	7.13
	Fe K $\beta$ 5	1	.03	1.7442	7.11

Nd L $\gamma$ 4	1	.1	1.7445	7.11
Gd L $\beta$ 2,15	1	20	1.7455	7.10
Tb L $\beta$ 3	1	6	1.7472	7.10
Sr K $\alpha$ 1,2	2	150	1.7533	14.14
Er Ln	1	1	1.7566	7.06
Fe K $\beta$ 1,3	1	20	1.7566	7.06

NOTES: Dy, Sm, Fe, Gd, and Sr can interfere with the Tm determination.  
Of these, only Sr can be excluded with a PHA.

YTTERBIUM

TABLE 70

ELEMENT: Yb  
 ANALYTICAL LINE: La1  
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 1.6719  
 BACKGROUND SETTING:  
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT % OXIDE</u>	<u>BACKGROUND SETTING</u>
Glass REE 2	4.26	±0.026

## INTERFERENCES:

<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
Ho L $\beta$ 1	1	50	1.6475	7.52
Tb L $\beta$ 5	1	.1	1.6510	7.51
Sm L $\gamma$ 3	1	1	1.6560	7.49
Eu L $\gamma$ 1	1	5	1.6574	7.48
Tb L $\beta$ 7	1	.1	1.6585	7.48
Ni K $\alpha$ 1,2	1	150	1.6592	7.47
Ho L $\beta$ 4	1	5	1.6595	7.47
Sm L $\gamma$ 2	1	1	1.6604	7.47
Y K $\alpha$ 1,2	2	150	1.6605	14.93
Tb L $\beta$ 10	1	.01	1.6673	7.44
YL La1	1	100	1.6719	7.41
Dy L $\beta$ 6	1	.1	1.6821	7.37
Dy L $\beta$ 3	1	6	1.6822	7.37
Yb La2	1	10	1.6829	7.37
Tb L $\beta$ 2,15	1	20	1.6830	7.37
Tm Ln	1	1	1.6963	7.31
Sm L $\gamma$ 6	1	.01	1.6967	7.31

NOTES: Of the interfering lines Ni K $\alpha$ 1,2; Y K $\alpha$ 1,2; DyL $\beta$ 3 and Tb<sub>L $\beta$ 2,15</sub> a PHA will only be useful to discriminate against Y. All other X-ray lines are of similar energy and impossible to exclude.



LUTETIUM

TABLE 71

ELEMENT: Lu  
 ANALYTICAL LINE: La1  
 CRYSTAL: LiF  
 SPECTROMETER SETTING: 1.6195  
 BACKGROUND SETTING:  
 ELEMENTAL SCANS:


<u>MINERAL</u>	<u>WT % OXIDE</u>	<u>BACKGROUND SETTING</u>
Glass REE 2	4.26	±0.026

## INTERFERENCES:


<u>LINE</u>	<u>N</u>	<u>I</u>	<u>LiF (200)</u>	<u>KeV</u>
Dy L $\beta$ 5	1	.1	1.5884	7.80
Yb L $_2$ -M $_2$	1	.01	1.5884	7.80
Eu L $\gamma$ 3	1	1	1.5903	7.80
Gd L $\gamma$ 1	1	5	1.5924	7.78
Eu L $\gamma$ 2	1	1	1.5961	7.77
Dy L $\beta$ 9	1	.01	1.5997	7.75
Er L $\beta$ 4	1	5	1.6007	7.74
Dy L $\beta$ 7	1	.1	1.6045	7.73
Sm L $\gamma$ 4	1	.1	1.6073	7.71
Dy L $\beta$ 10	1	.01	1.6074	7.71
Co K $\beta$ 5	1	.03	1.6089	7.71
Lu La1	1	100	1.6195	7.65
Co K $\beta$ 1,3	1	20	1.6208	7.65
Dy L $\beta$ 2,15	1	20	1.6237	7.63
Ho L $\beta$ 6	1	.1	1.6237	7.63

Eu L $\gamma$ 6	1	.01	1.6282	7.61
Lu La2	1	10	1.6303	7.60
Eu L $\gamma$ 8	1	.1	1.6346	7.58
Yb Ln	1	1	1.6356	7.58
Gd L $\gamma$ 5	1	.1	1.6412	7.55


NOTES: Chief interferences are from Co, Dy, Er and Gd. A PHA is of no help as an energy discriminator since all the lines are of similar X-ray energy.




La GLASS REE 3  
LiF 300 cps.  
La<sub>2</sub>O<sub>3</sub> 4.28 wt. %



Nd GLASS REE 2  
LiF 300 cps.  
Nd<sub>2</sub>O<sub>3</sub> 4.26 wt. %




Ce GLASS REE 3  
LiF 300 cps.  
Ce<sub>2</sub>O<sub>3</sub> 4.00 wt. %



68


FEET  
LEFT

REE IGWISI Hill Perovskite  
LiF 300 cps.






Pr GLASS REE 3  
LiF 300 cps.  
Pr<sub>2</sub>O<sub>3</sub> 4.45 wt. %




Sm GLASS REE 2  
LiF 300 cps.  
Sm<sub>2</sub>O<sub>3</sub> 4.26 wt. %







Eu GLASS REE 1  
LiF 300 cps.  
EuO 4.21 wt. %

Gd GLASS REE 1  
LiF 300 cps.  
Gd<sub>2</sub>O<sub>3</sub> 4.47 wt. %



Tb GLASS REE 1  
LiF 300 cps.  
Tb<sub>2</sub>O<sub>3</sub> 4.27 wt. %







Dy, Ho, Er	GLASS REE 4
LIF	300 cps.
Dy <sub>2</sub> O <sub>3</sub>	4.36 wt. %
Ho <sub>2</sub> O <sub>3</sub>	4.41 wt. %
Er <sub>2</sub> O <sub>3</sub>	4.36 wt. %




Tm GLASS REE 1  
LiF 300 cps.  
Tm<sub>2</sub>O<sub>3</sub> 4.36 wt. %



ORIGINAL PAGE IS  
OF POOR QUALITY





PLATINUM

TABLE 78

ELEMENT: Pt  
 CRYSTAL: PET  
 ANALYTICAL LINE:  $M\alpha$   
 SPECTROMETER SETTING: 2.7853  
 BACKGROUND SETTING:  $\pm 0.158$   
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT. %</u>	<u>BACKGROUND SETTING</u>
Pt	99.9	$\pm 0.158$


## INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	PET (002)	<u>KeV</u>
Pb	$M_3-N_1$	1	1	2.6273	2.17
Zr	$L\beta_6$	1	3	2.6301	2.17
Nb	$L\alpha_1$	1	100	2.6366	2.17
Mn	$K\beta_{1,3}$	3	20	2.6396	6.49
Nb	$L\alpha_2$	1	10	2.6401	2.16
P	$K\beta_1$	1	3	2.6733	2.14
Fe	$K\alpha_{1,2}$	3	150	2.6771	6.40
Pt	$M\beta$	1	50	2.6844	2.13
Zr	$L\beta_1$	1	45	2.6881	2.12
Au	$M\alpha_1$	1	100	2.6899	2.12
Th	$L\alpha_1$	3	100	2.7312	6.27
Y	$L\beta_3$	1	3	2.7559	2.07
Pt	$M\gamma-O_3$	1	.01	2.7576	2.07
Ni	$K\beta_{1,3}$	4	20	2.7639	8.26
Y	$L\beta_4$	1	3	2.7722	2.06
Pt	$M\alpha_1$	1	100	2.7853	2.05
Zr	$L\alpha_1$	1	100	2.7961	2.04
Y	$L\beta_6$	1	3	2.8070	2.03
P	$SK\alpha_4$	1	5	2.8137	2.03



P	SK $\alpha$ 3	1	5	2.8173	2.03
Mo	L1	1	3	2.8331	2.02
P	K $\alpha$ 1,2	1	150	2.8364	2.01
Ca	K $\beta$ 1,3	2	15	2.8463	4.01
Y	L $\beta$ 1	1	45	2.8613	2.00
Cr	K $\beta$ 1,3	3	18	2.8809	5.95
Au	M $\gamma$ -N $\gamma$	1	1	2.8829	1.98
Mn	K $\alpha$ 1,2	3	150	2.9061	5.89
Sr	L $\beta$ 3	1	3	2.9328	1.95

NOTES: With a material so widely used in experimental research it is difficult to decide which matrix will not present a problem and which will. The X-ray lines listed above should cover most of the geologic materials normally encountered. A PHA is useful in discriminating against some of the interfering lines — but not all of them. The following is a list of the lines where use of a PHA should be helpful: MnK $\beta$ 1,3, FeK $\alpha$ 1,2, ThL $\alpha$ 1, NiK $\beta$ 1,3, CaK $\beta$ 1,3, CrK $\beta$ 1,3, MnK $\alpha$ 1,2. For the following x-ray lines a PHA is of little or no help: NbL $\alpha$ 1, ZrL $\beta$ 1, AuM $\alpha$ 1, ZrL $\alpha$ 1, Pk $\alpha$ 1,2, YL $\beta$ 1.



ORIGINAL PAGE IS  
OF POOR QUALITY

THORIUM

TABLE 90

ELEMENT: Th  
 ANALYTICAL LINE: Th Ma1,2  
 CRYSTAL: PET  
 SPECTROMETER SETTING: 1.91 \_\_  
 BACKGROUND SETTING:  $\pm 0.020$

## ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT. % OXIDE</u>	<u>BACKGROUND SETTING</u>
GLASS W (PROBE SOCIETY)	0.79	$\pm 0.020$
Hibonite	0.53	$\pm 0.020$

## INTERFERENCES:

	<u>LINE</u>	<u>N</u>	<u>I</u>	<u>PET (002)</u>	<u>KeV</u>
	Gd La1	2	100	1.8855	6.06
	Ce Ly1	2	5	1.8873	6.05
	Bi M3-01	1	.5	1.8908	3.02
	Lu L $\beta$ 15	3	20	1.8952	9.04
	Gd La2	2	10	1.8957	6.02
	Tl M2-N4	1	1	1.8958	3.01
	Hf L $\beta$ 1	3	50	1.8987	9.02
	Rh L $\beta$ 2,15	1	25	1.9028	3.00
Th	Ma1	1	100	1.9060	3.00
	Pd L $\beta$ 1	1	42	1.9098	2.99
Th	Ma2	1	100	1.9120	2.99
	Ag La1	1	100	1.9135	2.98
	Pm L $\beta$ 1	2	50	1.9158	5.96
	Ag La2	1	10	1.9175	2.98
	Cr K $\beta$ 1,3	2	18	1.9206	5.95
	Ru Ly1	1	1	1.9263	2.96
	Ar Ka1,2	1	150	1.9312	2.96
	Cd Ln	1	1	1.9314	2.96
	Mn Ka1,2	2	150	1.9374	5.89

NOTES: There are several X-ray lines of similar energy as Th Ma and therefore not excluded with PHA. These include Bi, Tl, Rh, Pd, Ag, Ru, Ar and Cd. None of these however would typically be found in a Th matrix. Some REE might present a problem, but usually are of sufficiently low concentration that higher order reflections will not be seen.

URANIUM

TABLE 92

ELEMENT: U  
 ANALYTICAL LINE: UMa1,2  
 CRYSTAL: PET  
 SPECTROMETER SETTING: 1.80\_\_  
 BACKGROUND SETTING: +0.040  
 ELEMENTAL SCANS:

<u>MINERAL</u>	<u>WT. % OXIDE</u>	<u>BACKGROUND SETTING</u>
GLASS X (Probe Society)	0.76	±0.040

## INTERFERENCES:


<u>LINE</u>	<u>N</u>	<u>I</u>	<u>PET (002)</u>	<u>KeV</u>
Pd Lγ5	1	.1	1.7605	3.24
Th La1	4	100	1.7613	12.97
Ag Lβ3	1	11	1.7656	3.24
Bi M2-N4	1	5	1.7660	3.23
Dy La2	2	10	1.7686	6.46
Eu Lβ1	2	50	1.7690	6.64
Ag Lβ4	1	5	1.7826	3.20
Pb M1-N3	1	.1	1.7835	3.20
Fe Ka1,2	2	150	1.7847	6.40
Ar Kβ1,3	1	15	1.7899	3.19
Sb L1	1	7	1.7909	3.19
Bi M1-N2	1	1	1.7927	3.19
Ru Lγ2,3	1	.5	1.7953	3.18
Pd Lβ2,15	1	25	1.8004	3.17
U Ma1	1	100	1.8010	3.17
U Ma2	1	100	1.8074	3.16
Pb Lβ1	4	50	1.8109	12.61
Bi M3-04,5	1	1	1.8111	3.15



Ag L $\beta$ 1	1	42	1.8124	3.15
Zr K $\alpha$ 1,2	5	150	1.8132	15.74
Th M $\beta$	1	60	1.8152	3.15
Rh L $\gamma$ 1	1	10	1.8165	3.14
Tb L $\alpha$ 1	2	100	1.8208	6.27
Cd L $\alpha$ 1	1	100	1.8223	3.13
Pb M2-N4	1	5	1.8277	3.12
In Ln	1	7	1.8347	3.11
Tl M1-N3	1	1	1.8484	3.09
Pd L $\beta$ 6	1	1	1.8499	3.09

NOTES: All 1st order x-ray lines occurring within  $\pm Bkg$  are listed even though (1) background is only counted on the + side of peak and (2) many of the lines are not common to uranium minerals. Higher order lines for elements commonly occurring in uranium minerals are also given, i.e. Pb, Th, Fe, REE. X-ray lines from Th, REE, Fe, Zr and Pb L $\beta$ 1 can be excluded with a PHA; the remainder of x-ray lines cannot be discriminated against with a PHA. Most important of these lines is Th M $\beta$ .

Between oxides, carbonates, sulfates, phosphates, arsenates, vanadates, silicates, niobates, tantalates, titanates, and molybdates there are 100 uranium minerals listed in the "Handbook of Geochemistry" Vol. II-1, Springer-Verlag, 1969. Most commonly occurring uranium mineral is uraninite.



X-RAY OVERLAP AMONG TRANSITION ELEMENTS

TABLE 001

## X-RAY OVERLAP AMONG TRANSITION ELEMENTS

<u>ANALYTICAL LINE</u>	<u>INTERFERENCE</u>	<u>(Li F)</u> <u><math>\Delta\lambda</math></u>	<u><math>\Delta\text{KeV}</math></u>	<u>(PROBLEM)</u>
Ti $K\alpha_{1,2}$	Sc $K\beta$	0.0299	0.05	
V $K\alpha_{1,2}$	Ti $K\beta$	0.0091	0.02	(V in ilmenite)
Cr $K\alpha_{1,2}$	V $K\beta$	0.0066	0.02	
Mn $K\alpha_{1,2}$	Cr $K\beta$	0.0182	0.06	(Mn in chromite)
Fe $K\alpha_{1,2}$	Mn $K\beta$	0.0272	0.09	(Fe in rhodonite)
Co $K\alpha_{1,2}$	Fe $K\beta$	0.0337	0.14	(Co in kamacite)
Ni $K\alpha_{1,2}$	Co $K\beta$	0.0384	0.18	
Cu $K\alpha_{1,2}$	Ni $K\beta$	0.0417	0.22	

X-RAY OVERLAP AMONG RARE EARTH ELEMENTS

TABLE 002  
X-RAY OVERLAP AMONG RARE EARTH ELEMENTS

<u>Z</u>	<u>EL</u>	<u>LINE</u>	<u>LiF<math>\lambda</math></u>	<u>KeV</u>	<u>Z</u>	<u>EL</u>	<u>LINE</u>	<u>LiF<math>\lambda</math></u>	<u>KeV</u>
57	La	L $\alpha_1$	2.6657	4.65					
58	Ce	"	2.5615	4.84					
59	Pr	"	2.4630	5.03	57	La	L $\beta_1$	2.4589	5.04
60	Nd	"	2.3704	5.23	58	Ce	"	2.3561	5.26
61	Pm	"	2.2822	5.43					
62	Sm	"	2.1998	5.64					
63	Eu	"	2.1209	5.85					
64	Gd	"	2.0468	6.06					
65	Tb	"	1.9765	6.27	62	Sm	L $\beta_1$	1.9981	6.20
66	Dy	"	1.9088	6.49	63	Eu	"	1.9203	6.46
67	Ho	"	1.8450	6.72	64	Gd	"	1.8468	6.71
68	Er	"	1.7843	6.95	65	Tb	"	1.7768	6.98
69	Tm	"	1.7268	7.18	66	Dy	"	1.7106	7.25
70	Yb	"	1.6719	7.41					
71	Lu	"	1.6195	7.65					